



# ArtLiner™

— ORIGINAL —  
MICROPIGMENTS

ARTLINER - PMU colour number	ARTLINER colour name	Tone type (warm/cold/neutral)	Base colour	The best suitable Fitzpatrick skin type	ARTLINER MICROBLADING
101	Black Original	Neutral	Magenta/Blue/Violet	4,5,6	101-M
102	Dark Brown	Warm	Red/Orange/Yellow	3,4,5	102-M
103	Graphite	Neutral	Magenta/Blue/Violet	2,3	
104	Emerald Green	Neutral	Magenta/Blue/Violet	1,2	
105	Smaragd Green	Neutral	Magenta/Blue/Violet	1,2,3	
106	Aquamarine	Neutral	Magenta/Blue/Violet	2,3,4	
107	Denim Blue	Neutral	Magenta/Blue/Violet	2,3	
108	Pacific Blue	Neutral	Magenta/Blue/Violet	2,3	
109	Navy Blue	Neutral/Cold	Magenta/Blue/Violet	2,3	
110	Aubergine	Cold	Magenta/Blue/Violet	3,4	
201	Black Brown	Neutral	Red/Orange/Yellow	4,5,6	201-M
202	Mokka	Warm	Red/Orange/Yellow	3,4	202-M
203	Medium Brown	Warm	Red/Orange/Yellow	2,3	203-M
204	Tropical Brown	Warm	Red/Orange/Yellow	3,4,5	204-M
205	Taupe	Cold	Magenta/Blue/Violet	2,3	205-M
206	Mystic Brown	Warm	Red/Orange/Yellow	1,2,3	206-M
207	Monroe	Warm	Red/Orange/Yellow	1,2	207-M
208	Hazelnut	Warm	Red/Orange/Yellow	1,2	208-M
209	Golden Sand	Warm	Red/Orange/Yellow	1,2	209-M
210	Blond	Warm	Red/Orange/Yellow	1,2	210-M
301	Pink Rose	Cold	Magenta/Blue/Violet	For the lipcolors it's difficult to select any skintype. The lipcolor depends of the natural lipcolor. Some people have a blue undertone, other one have a rose undertone or a purple. The best is to mix 2 or 3 colors for the lips.	
302	Amaranth	Warm	Magenta/Blue/Violet		
303	Burgundy	Cold	Magenta/Blue/Violet		
304	Dark Rose	Warm	Red/Orange/Yellow		
305	Coral	Warm	Red/Orange/Yellow		
306	Orange Red	Warm	Red/Orange/Yellow		
307	Simply Red	Warm	Red/Orange/Yellow		
308	Wine Red	Warm	Red/Orange/Yellow		
309	Red Brown	Neutral	Red/Orange/Yellow		
310	Island Brown	Warm	Red/Orange/Yellow		
401	Pale	Neutral	Red/Orange/Yellow	For corrective pigments it's also not possible to select a skin type because they are also used for medical pigmentation at breast. So even a person with skin type 1 can sometimes need for exemple n°406 (Tan) for breast pigmentation.	
402	Caucasian Skin	Warm	Red/Orange/Yellow		
403	Skin Light	Warm	Red/Orange/Yellow		
404	Skin Medium	Warm	Red/Orange/Yellow		
405	Skin Dark	Neutral/Warm	Red/Orange/Yellow		
406	Tan	Warm	Red/Orange/Yellow		
407	Skin Rose	Cold	Magenta/Blue/Violet		
408	Vanilla	Warm	Red/Orange/Yellow		
409	Ivory	Warm	Red/Orange/Yellow		
410	White	Neutral	Magenta/Blue/Violet		

[www.art-liner.com](http://www.art-liner.com)